



StabilCoat™ and 5X StabilCoat Immunoassay Stabilizer

Product insert

Product number:

SC01-0050 (50 mL)
SC01-1000 (1000 mL)
SC01-2000 (2000 mL)

SC05-0100 (100 mL)
SC05-1000 (1000 mL)
SC05-2000 (2000 mL)

Product use:

StabilCoat™ Immunoassay Stabilizer, and the 5X concentrate, effectively preserves the conformation and activity of dried proteins while simultaneously blocking the surface to reduce non-specific binding. The product can be used to stabilize antibodies, antigens or enzymes on an assortment of immunoassay components: polystyrene plates, tubes, glass, membranes, and filter paper.

The non-specific binding of serum proteins, conjugates and other assay components to solid surfaces in immunoassay applications requires that a blocking step be incorporated into the protocol to maximize assay sensitivity. Additionally, the manufacturing of immunoassays test kits requires long term protein stability to increase the manufacturing scale and extended shelf lives for the end user. The StabilCoat formulation allows for these requirements to be accomplished in a single step while easily being substituted for the blocking solution in an existing assay protocol.

Product stability, storage and specifications:

Product stability	StabilCoat: 3 years from date of manufacture 5X StabilCoat: 2 years from date of manufacture
Storage	Product should be stored at 2-8° C or room temperature
Specification	Protein: Bovine Product Buffer: PBS pH: 7.0-7.4 when at working solution Preservative: None
Notes	Please note that Surmodics™ Stabilization Products are shipped to customers at ambient temperature. Extensive stability studies have shown that prolonged storage at ambient temperature will not affect the product quality or efficacy.

Recommendations for use:

Aseptically pour off the desired volumes of material needed for the application and allow the product to equilibrate at room temperature prior to use.

The following are general guidelines only.

To Prepare a Dilution of 5X StabilCoat Stabilizer:

- 1) Add one part 5X StabilCoat to four parts deionized water.
- 2) Mix the solution without foaming to obtain a homogenous mixture.

StabilCoat™ and 5X StabilCoat Immunoassay Stabilizer

Product insert

To Stabilize Adsorbed or Immobilized Proteins on Microtiter Plates/Strips:

- 1) Immobilize or adsorb the primary protein (antibody or antigen) to the plate surface according to an optimized method.
- 2) Wash the plate to remove excess or weakly bound protein.
- 3) Add StabilCoat to allow interaction with the entire protein-coated surface.
 - For example: if the total assay volume is 150 µL, add at least 200 µL of StabilCoat to each well.
- 4) Incubate up to 60 minutes at room temperature. Reduced incubation times can be achieved in certain applications.
- 5) Aspirate the excess StabilCoat from the microwells.
- 6) Dry components for long-term storage using a method that is optimized by the end user. Recommended methods are as follows:
 - Dry in a humidity controlled chamber (less than 15% humidity) until dry (4 to 24 hours)
 - Dry 2-4 hours in a sealed desiccated container.
 - Dry plates at 30-40°C in a vacuum oven for 4 hours.
- 7) Package the stabilized coated plate in a heat-sealed moisture proof pouch with a desiccant.

Additional considerations:

A preservative may be added to StabilCoat if required by the assay protocol.

For technical assistance, email ivdtechsupport@surmodics.com

Related products:

Blockers/Stabilizers:
StabilGuard™ Immunoassay Stabilizer (BSA-Free) (SG01)
StabilGuard™ Choice Stabilizer(BSA-Free) (SG02)
StabilCoat™ Plus Stabilizer (SC02)
StabilBlock™ Immunoassay Stabilizer (ST01)
In-Solution Protein Stabilizers & Diluents:
StabilZyme™ Protein Free Stabilizer (SZPF)
StabilZyme™ Protein Free AP Stabilizer (SAPF)
StabilZyme™ AP Conjugate Stabilizer (SA01)
StabilZyme™ HRP Conjugate Stabilizer (SZ02)
StabilZyme™ SELECT Stabilizer (SZ03)
StabilZyme™ NOBLE Stabilizer (SZ04)
Surmodics™ Assay Diluent (Protein-Free) (SM01)

StabilCoat™ and 5X StabilCoat Immunoassay Stabilizer

Product insert

Substrates:
BioFX™ TMB Supersensitive One Component HRP Microwell Substrate (TMBS)
BioFX™ TMB One Component HRP Microwell Substrate (TMBW)
BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)
BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)
BioFX™ TMB Extended Range HRP Microwell Substrate (TMBX)
BioFX™ TMB Enhanced HRP Membrane Substrate (ESPM)
BioFX™ TMB One Component HRP Membrane Substrate (TMBM)
BioFX™ pNPP Yellow One Component Microwell Substrate with Stabilizing Pellets (PNPS)
BioFX™ ABTS One Component HRP Microwell Substrate (ABTS)

Products distributed by



Part of  **BBI Solutions**

DIARECT GmbH
Bötzingen Str. 29 B
79111 Freiburg
Germany

Tel. +49 (0) 761 47979-0
Fax +49 (0) 761 47979 -29

www.bbisolutions.com
orders-dia@bbisolutions.com
info-dia@bbisolutions.com

SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD and STABILBLOCK are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for further manufacture and/or to be used as a component with in-vitro diagnostics immunoassays and are not intended for use in human or therapeutics purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made according to DIARECT's Terms and Conditions of Sale.

Aug-2020

Product Code: SC01, SC05
Doc # 7606, Revision 03

Page 3 of 3