



StabilZyme™ Protein-Free Stabilizer

Product insert

Product number:

SZPF-0050 (50 mL)
SZPF-1000 (1000 mL)
SZPF-2000 (2000 mL)

Product use:

StabilZyme™ conjugate/protein stabilizers protect the entire conjugate by preventing the loss of catalytic activity and maintaining the structural integrity of the protein in solution. These stabilizers allow for the storage of conjugated proteins at lower use concentrations, extended shelf life and increased signal-to-noise ratios for improved assay performance.

The StabilZyme Protein-Free Stabilizer is optimized to maintain the conformation of antibodies/antigens conjugated to horseradish peroxidase (HRP), antibody-coated particles and unmodified proteins in solution. The key benefits of the protein-free formulation are the elimination of protein interference and cross reactivity with the purpose of reducing erroneous results within immunoassay applications.

Adding the conjugate directly into the StabilZyme Protein-Free Stabilizer allows for stable storage of the conjugate at a working concentration and eliminates the need for subsequent dilutions.

Product stability, storage and specifications:

Product stability	Stable for 2.5 years from date of manufacture
Storage	Product should be stored at 2-8°C or at room temperature
Specification	Protein: None Product Buffer: MOPS 6.3-7.3 Preservative: 0.02% methylisothiazolone and 0.02% bromonitrodioxane
Notes	Please note that Surmodics™ Stabilization Products are shipped to customers at ambient temperature. Extensive stability studies have shown that prolonged storage at ambient temperature will not affect the product quality or efficacy.

Recommendations for use:

Aseptically pour off desired volumes of material needed for the application and allow the product to equilibrate at room temperature prior to use.

The following are general guidelines only.

- 1) Use StabilZyme Protein-Free Stabilizer at 100% concentration for optimal stability and performance.
 - If necessary, adjust the pH of the solution to align with the properties of the protein for optimal performance.
- 2) Dilute the HRP conjugate/protein/antibody coated particles to a working concentration in the StabilZyme Protein Free Stabilizer solution.
- 3) Use the diluted conjugate solution according to the standard assay protocol.
- 4) For optimal performance, store the conjugate solution at 2-8°C and protect from direct exposure to light.

For technical assistance, email ivdtechsupport@surmodics.com

StabilZyme™ Protein-Free Stabilizer

Product insert

Related products:

In-Solution Protein Stabilizers & Diluents:
StabilZyme™ HRP Conjugate Stabilizer (SZ02)
StabilZyme™ SELECT Stabilizer (SZ03)
StabilZyme™ NOBLE Stabilizer (SZ04)
Surmodics™ Assay Diluent (Protein-Free) (SM01)
MatrixGuard™ Assay Diluent (SM02)
Blockers/Stabilizers:
StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)
StabilCoat™ Immunoassay Stabilizer (SC01)
StabilBlock™ Immunoassay Stabilizer (ST01)
Substrates:
BioFX™ TMB Super Sensitive One Component HRP Microwell Substrate (TMBS)
BioFX™ TMB One Component HRP Microwell Substrate (TMBW)
BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)
BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)
BioFX™ TMB Extended Range HRP Microwell Substrate (TMBX)
BioFX™ TMB Enhanced HRP Membrane Substrate (ESPM)
BioFX™ TMB One Component HRP Membrane Substrate (TMBM)
BioFX™ ABTS One Component HRP Microwell Substrate (ABTS)

Products distributed by



DIARECT GmbH
Bötzingen Str. 29 B
79111 Freiburg
Germany

Tel. +49 (0) 761 47979-0
Fax +49 (0) 761 47979 -29

www.bbisolutions.com
orders-dia@bbisolutions.com
info-dia@bbisolutions.com

SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD and STABILBLOCK are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for further manufacture and/or to be used as a component with in-vitro diagnostics immunoassays and are not intended for use in human or therapeutics purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made according to DIARECT's Terms and Conditions of Sale.