

Product number:

TMBM-0100-01 (100 mL)
 TMBM-1000-01 (1000 mL)
 Custom Sizes Available

Product use:

BioFX™ TMB One Component HRP Membrane Substrate is a ready-to-use solution containing 3,3',5,5'-tetramethylbenzidine in a mildly acidic buffer. When this substrate system is reacted with peroxidase, an insoluble permanent dark blue reaction product is obtained.

Product stability, storage and specifications:

Product stability	Stable for 2 years from date of manufacture when stored at 15-25°C. Stable for 4 years when stored at 2-8°C.
Storage	This substrate is light sensitive and should be protected from direct ultraviolet sources
Specification	Unreacted substrate should be faint yellow to pink in appearance. Over time, some flocculent material may appear; however, this does not affect product performance

Recommendations for use:

The product should be allowed to equilibrate to room temperature (25°C) prior to use.

- 1) Cover the membrane surface with an ample amount of substrate and incubate for 5-20 minutes.
 - The substrate system will react with sites on the membrane containing peroxidase producing an insoluble permanent dark blue reaction product.
 - The reaction should be monitored and read before the background color becomes too intense which may result in diminished contrast between the positive and background staining.
- 2) Stop the reaction by rinsing the membrane with reagent quality water.


Additional considerations:

Dilution of the substrate is not recommended. It is recommended that the antibodies or conjugates be diluted to reduce the intensity of the reaction.

For technical assistance, email ivdtechsupport@surmodics.com

Related products:

In-Solution Protein Stabilizers & Diluents:
StabilZyme™ HRP Conjugate Stabilizer (SZ02)
StabilZyme™ SELECT Stabilizer (SZ03)
StabilZyme™ NOBLE Stabilizer (SZ04)
StabilZyme™ Protein Free Stabilizer (SZPF)
Surmodics™ Assay Diluent (Protein-Free) (SM01)
MatrixGuard™ Diluent (SM02)

 **SURMODICS** **BioFX™ TMB One Component HRP Membrane Substrate**

Product insert

Blockers/Stabilizers:
StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)
StabilCoat™ Immunoassay Stabilizer (SC01)
StabilBlock™ Immunoassay Stabilizer (ST01)
StabilBlot™ Protein-Free Blocker (SB01)
StabilBlot™ Milk Blocker – 10X (SB04)
StabilBlot™ Casein Blocker – 5X (SB05)
Substrates:
BioFX™ Chemiluminescent Ultra-Sensitive HRP Microwell Substrate (CHMI)
BioFX™ Chemiluminescent Sensitive Plus HRP Microwell and/or Membrane Substrate (LERI)
BioFX™ Chemiluminescent Super Sensitive HRP Microwell and/or Membrane Substrate (LUMI)
BioFX™ TMB Super Sensitive One Component HRP Microwell Substrate (TMBS)
BioFX™ TMB One Component HRP Microwell Substrate (TMBW)
BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)
BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)
BioFX™ TMB Extended Range One Component HRP Microwell Substrate (TMBX)
BioFX™ TMB Enhanced HRP Membrane Substrate (ESPM)

Products distributed by



DIARECT GmbH
 Bötzingen Str. 29 B
 79111 Freiburg
 Germany

Tel. +49 (0) 761 47979-0
 Fax +49 (0) 761 47979 -29

www.bbisolutions.com
 orders-dia@bbisolutions.com
 info-dia@bbisolutions.com

SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD, STABILBLOCK and STABILBLOT are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for further manufacture and/or to be used as a component with in-vitro diagnostics immunoassays and are not intended for use in human or therapeutics purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made according to DIARECT's Terms and Conditions of Sale.

19-May-2020

Product Code: TMBM
 Doc # 7775, Revision 01

Page 2