

<b>Originating Department</b>	QC
<b>Approval Departments</b>	QC & Validation
<b>Effective Date</b>	Refer to Q-Pulse

## 1.0 PRODUCT DETAILS


1.1 **Enzyme Name:** Catalase

1.2 **Systematic Name:** Hydrogen-peroxide:hydrogen-peroxide oxidoreductase

1.3 **E.C. Number:** 1.11.1.6

1.4 **Source:** Bovine liver & *Aspergillus niger*

## 2.0 ASSAY PRINCIPLE



The decrease in absorbance at 240nm is proportional to the catalase activity.

## 3.0 UNIT DEFINITION

That amount of enzyme causing the decomposition of one micromole of hydrogen peroxide per minute at 25°C and pH 7.0.

## 4.0 EQUIPMENT REQUIRED

Double beam UV/vis spectrophotometer with chart recorder  
Water bath set to achieve a reaction temperature of 25°C (± 0.1°C)  
Thermometer  
Silica cuvettes  
Test tubes  
Manual pipettes and tips

## 5.0 REAGENTS REQUIRED

When using hazardous chemicals, handle in accordance with COSHH Regulations.

## Reagent details

Chemical / Reagent	Supplier	Product No.	F.W.
6M Hydrochloric acid	Fisher Scientific	72033-1L	N/A
Di-sodium hydrogen phosphate dihydrate	VWR	28029.260	177.99
Potassium dihydrogen phosphate	VWR	26936.293	136.09
Hydrogen peroxide solution	Sigma	1.07209.0250	N/A

## 6.0 PREPARATION OF REAGENTS

### 6.1 6M Hydrochloric acid

Use as required and refer to the manufacturer's expiry date.

### 6.2 0.05M Sodium potassium phosphate pH7.0

Dissolve 5.87g of di-sodium hydrogen phosphate dihydrate and 2.28g of potassium dihydrogen phosphate in approximately 900ml of water. Adjust to pH 7.0 with 6M Hydrochloric acid and make up to 1 litre with water.

Stable for 2 weeks at 2 to 8°C

### 6.3 Working substrate

Add approximately 0.64ml<sup>1</sup> of 30% hydrogen peroxide to 100ml of 0.05M Sodium potassium phosphate pH 7.0.

Pipette 2ml of 0.05M Sodium potassium phosphate pH 7.0 and 1ml of working substrate into a test tube. Mix thoroughly and transfer to a silica cuvette. Measure the absorbance at 240nm versus 0.05M Sodium potassium phosphate pH 7.0. If necessary, add more 0.05M Sodium potassium phosphate pH 7.0 or 30% hydrogen peroxide to the working substrate to adjust the concentration until an absorbance of 0.85 (±0.02) is obtained.

Store in a dark bottle. Stable for 1 month when stored at 2 to 8°C.

### 6.4 Enzyme solution

Freeze-dried powders:

Accurately weigh at least 10mg of freeze-dried powder into new glass vials, each test sample to be weighed in triplicate. Dissolve each to a concentration of 5mg/ml in 0.05M Sodium potassium phosphate pH 7.0. Immediately prior to assay, dilute to approximately 0.6U/ml in 0.05M Sodium potassium phosphate pH 7.0, the final 4.5ml serial dilution buffer having been **pre-equilibrated at 25°C.**<sup>2</sup>

<sup>1</sup> 0.64ml is a guide as it can vary between different lots and over time.

<sup>2</sup> Pre-equilibration is required since the volume of diluted enzyme (2ml) is a significant proportion of the reaction mixture volume (3ml) and hence will influence the reaction temperature.

Liquid preparations:

Immediately prior to assay, dilute to approximately 0.6U/ml in 0.05M Sodium potassium phosphate pH 7.0, the final 4.5ml serial dilution buffer having been **pre-equilibrated at 25°C.**<sup>2</sup>

Process samples:

Dilute to approximately 0.6U/ml in 0.05M Sodium potassium phosphate pH 7.0, ensuring the concentration is within the range of 0.0297 U/ml to 0.914 U/ml (equivalent to reaction rates ( $\Delta A_{240}/\text{min}$ ) of 0.0086 to 0.0266).<sup>3</sup>

## 7.0 PROCEDURE

Allow 0.05M Sodium potassium phosphate pH 7.0 and working substrate to equilibrate at 25°C for at least 5 minutes in the water bath before use.

Temperature = 25°C

Wavelength = 240nm

Light path = 10mm

Into disposable test tubes pipette the following:

	<b>Test</b>	<b>Reference</b>
Working substrate	1.00ml	0.00ml
0.05M Sodium potassium phosphate pH 7.0	0.00ml	3.00ml

Allow the solutions to equilibrate to 25°C then add:

Enzyme solution, diluted to ~0.6U/ml:	<u>2.00ml</u>	<u>0.00ml</u>
Total reaction mix volume ( $V_t$ ):	3.00ml	3.00ml

Transfer the solutions to silica cuvettes then record the decrease in absorbance at 240nm, reading the test solution versus the reference solution for approximately 3 minutes. Measure the change in absorbance per minute ( $\Delta A_{240}/\text{min}$ ) over the linear portion of the curve and use this value in the calculation.

---

<sup>3</sup> Taken from Analytical Test Method Validation ATMV-013

## 8.0 CALCULATION

$$8.1 \text{ Volume activity (U/ml)} = \frac{\Delta A_{240}/\text{min} \times V_t \times \text{dilution factor}}{V_s \times \epsilon}$$

Where:  $V_t$  = final volume of the reaction mix (ml) = 3.00  
 $V_s$  = sample volume (ml) = 2.00  
 $\epsilon$  = micromolar extinction coefficient for peroxide ( $\text{cm}^2/\mu\text{mole}$ ) = 0.0436

$$\text{Volume activity (U/ml)} = \Delta A_{240}/\text{min} \times 34.4 \times \text{dilution factor}$$

$$8.2. \text{ For freeze-dried samples: Weight activity (U/mg material)} = \frac{\text{U/ml}}{\text{mg/ml}}$$

$$\text{Specific activity (U/mg protein)} = \frac{\text{U/mg material}}{\text{mg protein/mg material}}$$

8.3 For liquid preparations:

$$\text{Specific activity (U/mg protein)} = \frac{\text{U/ml}}{\text{mg protein/ml}}$$

## 9.0 PROTEIN DETERMINATION

Protein is determined by the method of Lowry et al in accordance with Analytical Procedure AP62<sup>2</sup>.

## 10.0 $A_{280}^{1\%}$ DETERMINATION

This is determined in accordance with Analytical Procedure AP63.

## 11.0 ASSOCIATED DOCUMENTS

AP62	Lowry Protein Determination
AP63	Spectrophotometric Measurements
ATMV-013	Validation of Catalase (Kinetic Method)

## 12.0 REFERENCES

1. Valle, B. L. and Hoch F. L. (1955) *Proc. Nat. Acad. Sci.* **41**, 327
2. Lowry, O.H., Rosebrough, N.J., Farr, A.L., & Randall, R.J. (1951) *J. Biol. Chem.* **193**, 265

**REVISION HISTORY**

Document version number	Section number	Summary of Changes
09	Global	Header Changed to reflect current practice. Approval date removed and Effective Date changed to 'refer to Q-Pulse'. Removed document history summary for 08 to provide space for further changes.
	5.0	Reagent details- changed suppliers to Sigma and corrected product numbers from hydrochloric acid and hydrogen peroxide.
	7.0	Changed Enzyme solution, diluted to ~0.6U/ml from ~0.2U/ml in line with enzyme solution dilution in section 6.4.