

Product Analysis

EM Rabbit anti Goat IgG (H&L): 15nm Gold

Product Information

Product Code:	EM.RAG15
Batch Number:	15793
Protein Absorption:	Not Absorbed
Protein Concentration:	8.7µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<1.5ng	Passed
Optical Density:	4.0	Passed
Mean Diameter:	14.1nm	Passed
Coefficient of Variation:	6.0%	Passed
Percentage of Single Particles:	97%	Passed

For in vitro research only. Not for diagnostic use