

SRP54 humAb IgG

Antibody Specification

Product Number: 36500

Host: GammaPrim[®] mouse

Clonality: monoclonal

Isotype: IgG 1 / kappa

Description:

Human chimeric monoclonal antibody raised against recombinant SRP54.

For in vitro research and manufacturing use only.

Immunogen/Epitope:

Human SRP54.

Immunological function:

ELISA and Membrane-based immunoassay with DIARECT's SRP54 antigen.

Origin:

Produced by heterohybridoma cell line (constructed by fusing peripheral B lymphocytes from a GammaPrim[®] mouse immunized with SRP54 to a hetero-hybridoma fusion cell line).

Concentration for assay use:

Concentrations for use as standard, calibrator and/or positive control reagent in immunodiagnostic assays have to be individually determined by user.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral pH and 10% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2019 - 2021 DIARECT GmbH