

Product Analysis

EM Goat anti-Mouse IgG (H&L): 10nm Gold

Product Information

Product Code:	EM.GMHL10
Batch Number:	*****
Protein Absorption:	Minimum cross reaction to Human, Bovine and Horse serum proteins
Protein Concentration:	*** µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<***pg after silver enhancement	Passed
Optical Density:	*****	Passed
Mean Diameter:	*nm	Passed
Coefficient of Variation:	***%	Passed
Percentage of Single Particles:	**%	Passed

For in vitro research only. Not for diagnostic use