

## **Product Analysis**

### **EM Protein A: 15nm Gold**

#### **Product Information**

Product Code:	EM.PAG15
Batch Number:	*****
Protein Absorption:	Not Applicable
Protein Concentration:	***µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

#### **Quality Control**

Protein Detected:	<***pg	<b>Passed</b>
Optical Density:	*****	<b>Passed</b>
Mean Diameter:	*nm	<b>Passed</b>
Coefficient of Variation:	***%	<b>Passed</b>
Percentage of Single Particles:	**%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**