

Product Analysis

EM Streptavidin: 10nm Gold

Product Information

Product Code:	EM.STP10
Batch Number:	*****
Protein Absorption:	Not Applicable
Protein Concentration:	*** µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<***pg after silver enhancement	Passed
Optical Density:	*****	Passed
Mean Diameter:	*nm	Passed
Coefficient of Variation:	***%	Passed
Percentage of Single Particles:	**%	Passed

A pdf protocol for the use of Streptavidin conjugates is available on request

For in vitro research only. Not for diagnostic use

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483