

Thyroid Peroxidase (TPO) biotinylated

Antigen Specification

Product Number: 20100

Description:

Human thyroid peroxidase (TPO); identical with the antigen formerly called “microsomal antigen”. Biotinylated recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human auto-antibodies. Auto-antibodies to TPO recognize conformation-dependent epitopes.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Truncated cDNA coding for the extracellular part of human TPO. Transmembrane and intracytoplasmic parts of the molecule are removed and replaced by a hexa-histidine purification tag at the C-terminus. After cleavage of the secretory leader pre-sequence the mature recombinant TPO has a length of 834 amino acids.

Biochemical tests:

Electrophoresis (purity > 80%)

Calculated molecular weight:

93 kDa (for protein component excluding glycosylation).

Calculated isoelectric point:

pH 6.2

Immunological tests/Functionality:

Functional Streptavidin-based ELISA test (analysis of positive/negative samples).

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: ionic strength around 150 mM, neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided. Use of reducing sulfhydryl reagents and denaturing detergents will result in the destruction of conformation epitopes.

Copyright 2017 - 2021 DIARECT GmbH