

## Ara h 9.0101

### Allergen Specification

**Product Number:** 52000

**Description:**

*Arachis hypogaea* Ara h 9.0101.  
Recombinant allergen for in vitro  
research and manufacturing use only.

**Immunological function:**

Binds IgE-type human antibodies.

**Origin:**

Recombinant. Expressed by  
recombinant baculovirus (*Autographa  
californica* multiple nuclear  
polyhedrosis virus; AcMNPV) infection  
of *Spodoptera frugiperda* Sf9 insect  
cells.

**Expression construct:**

Full-length cDNA coding for Ara h  
9.0101 fused to a hexa-histidine  
purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western  
blot with patient sample.

**Calculated molecular weight:**

10 kDa

**Calculated isoelectric point:**

pH 9.45

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte  
Lysate) chromogenic endotoxin assay.

**Immunological tests/Functionality:**

Immunodot analyses with  
positive/negative samples.

**Recommended buffer/storage and  
handling conditions:**

Recommendations for storage buffer:  
neutral to slightly alkaline pH and 20 %  
glycerol as cryoprotective agent.  
Storage conditions: -70°C or below.  
Repeated freeze/thaw cycles should  
be avoided.

Copyright 2017 - 2021 DIARECT GmbH