

## Phl p 12.0101

### Allergen Specification

**Product Number:** 52400

**Description:**

*Phleum pratense* Phl p 12.0101.  
Recombinant allergen for in vitro  
research and manufacturing use only.

**Immunological function:**

Binds IgE-type human antibodies.

**Origin:**

Recombinant. Expressed by  
recombinant baculovirus (*Autographa  
californica* multiple nuclear  
polyhedrosis virus; AcMNPV) infection  
of *Spodoptera frugiperda* Sf9 insect  
cells.

**Expression construct:**

Full-length cDNA coding for Phl p  
12.0101 fused to a deca-histidine  
purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western  
blot with patient sample.

**Calculated molecular weight:**

16 kDa

**Calculated isoelectric point:**

pH 6.06

**Molar extinction coefficient:**

17085;  $A_{280}$  (1 mg/mL) = 1.095

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte  
Lysate) chromogenic endotoxin assay.

**Immunological tests/Functionality:**

Immunodot analyses with  
positive/negative samples.

**Recommended buffer/storage and  
handling conditions:**

Recommendations for storage buffer:  
neutral to slightly alkaline pH and 20%  
glycerol as cryoprotective agent.  
Storage conditions: -70°C or below.  
Repeated freeze/thaw cycles should  
be avoided.

Copyright 2017 - 2021 DIARECT GmbH