

<b>PROCESS:</b> Product Specifications	<b>ISSUE DATE:</b> 01.10.2016	<b>DOC NO:</b> 161905BBI
<b>TITLE:</b> Trypsin/Chymotrypsin 6:1	<b>REVISION:</b> 5	<b>PAGE:</b> 1 of 1
<b>AUTHORISED:</b> Product Manager		<b>APPROVED:</b> QA Manager

Product Code: 161905BBI EC no.:3.4.21.4 EDQM Certificate no.: R1-CEP 2001-209-Rev 03

Aspect	Description	Units	Limits	Method
<b>Form</b>	Off-white, freeze-dried powder	N/A	N/A	N/A
<b>Preparation</b>	Bovine pancreas	N/A	N/A	N/A
<b>Activity</b>	Trypsin Chymotrypsin	U/mg material U/mg material	2200-2760 365-460	BP6.2.2-WI20 BP6.2.2-WI12
<b>Retest Date</b>	Re-assay	Months	36	N/A
<b>Storage</b>	Storage Temperature	°C	-20	N/A
<b>Transit</b>	Storage in Transit	°C	Ambient	N/A
<b>Packaging</b>	Dependent on volume / weight	N/A	N/A	N/A
<b>Labels</b>	<u>Large Label</u> <ul style="list-style-type: none"> <li>• Name of Product</li> <li>• Product Code</li> <li>• Lot Number</li> <li>• Net Weight</li> <li>• Manufacturing Date</li> <li>• Retest Date</li> <li>• Storage Temperature</li> <li>• Storage in Transit</li> </ul> <u>Small Label</u> <ul style="list-style-type: none"> <li>• Name of Product</li> <li>• Product Code</li> <li>• Lot No.</li> <li>• Net Weight</li> <li>• Manufacture Date</li> <li>• Retest Date</li> <li>• Storage Temperature</li> <li>• Storage in Transit</li> </ul>	N/A  °C °C	N/A  -20 Ambient	N/A
<b>Samples</b>	Retention	1 g	3	N/A
	Quality Control	1 g	1	
<b>Special Instructions</b>	None			