

Product:	ALKALINE PHOSPHATASE
Product code:	ALPI11G
E.C. number:	3.1.3.1
CAS number:	9001-78-9
EINECS number:	232-631-4
Systematic name:	Orthophosphoric-monoester phosphohydrolase (alkaline optimum).
Alternative name:	Alkaline phosphomonoesterase; Phosphomonoesterase; Glycerophosphatase.
Source:	Bovine intestinal mucosa
Form:	A clear, almost colourless solution in 40% Glycerol containing 6mM Tris/HCl, 6mM Magnesium Chloride and 0.12mM Zinc Chloride, pH approximately 7.6
Storage conditions:	Store at 2°C to 8°C. <i>DO NOT FREEZE.</i>
Unit definition:	<i>Glycine:</i> That amount of enzyme causing the hydrolysis of one micromole of <i>p</i> -nitrophenyl phosphate per minute at 25°C and pH 9.6 (glycine buffer). <i>DEA:</i> That amount of enzyme causing the hydrolysis of one micromole of <i>p</i> -nitrophenyl phosphate per minute at 37°C and pH 9.8 (diethanolamine buffer).
Specific activity:	Not less than 1800 Glycine U/mg protein , (Equivalent to approximately 5400 DEA U/mg protein).
Bovine IgG:	Not detectable
Chromatogram analysis:	Greater than 90% pure by molecular exclusion chromatography
Typical properties	
Protein Concentration:	10-20 mg/ml (Determined by Biuret procedure)
Available amino groups:	8 - 13 moles amino groups per mole enzyme
Carbohydrate content:	4.0 – 6.5%