

<b>Product:</b>	<b>TRYPSIN INHIBITOR</b>
<b>Product code:</b>	<b>SOY1</b>
<b>CAS number:</b>	9035-81-8
<b>EINECS number:</b>	232-906-9
<b>Source:</b>	Soybean
<b>Form:</b>	A freeze-dried material
<b>Storage conditions:</b>	Store desiccated at -15°C or below. Allow to come to room temperature before opening. Before returning to storage, re-desiccate under vacuum over silica gel for a minimum of four hours. Re-seal before returning to -15°C or below.
<b>Unit definition:</b>	That amount of inhibitor that will inhibit one unit of Trypsin activity (BAEE Units)
<b>Activity:</b>	Not less than <b>7,000 BAEE Units of inhibition/mg material</b>
<b>Solubility:</b>	Dissolves readily at 5mg/ml in 0.067M Potassium phosphate buffer, pH 7.6 to give a clear solution.