

MANUAL: Specification	ISSUE DATE: 29.11.2021	DOC NO: 121705BBI Q-pulse Ref. No.: CT-PRDS16
TITLE: Deoxyribonuclease 20/40 % A/S ppt	REVISION: 2	PAGE: 1 of 1
AUTHORISED: Laboratory Supervisor		APPROVED: SQR Leader

Product Code: 121705BBI

EC no.: N/A

Aspect	Description	Units	Limits	Method
Form	Tan-white 20/40% A/S precipitate.	N/A	N/A	N/A
Preparation	Bovine pancreas	N/A	N/A	N/A
Activity	DNase Activity (Kunitz) 10mg/ml	u/mg material	Record	BP6.2.2-WI13
Additional Data	E ₂₈₀ 1% (10mg/ml)	Absorbance	Record	BP6.2.2-WI138
	Specific Activity	U/A280	Record	N/A
	Total Aerobic Microbial Count	cfu/g	Record	BP6.2.3-WI49
Storage	Storage Temperature	°C	-20	N/A
Transit	Storage in Transit	°C	Ambient	N/A
Packaging	<ul style="list-style-type: none"> 15g Sample to be submitted with each order. Double polyethylene bag in a drum. 	N/A	N/A	N/A
Retest Date	Re-assay	Months	24	N/A
Labels	<u>Large Labels for bulk material</u> <ul style="list-style-type: none"> Name of Product Product Code Lot No. Net Weight Manufacture Date Retest Date Storage Temperature Storage in Transit 	N/A °C °C	N/A -20 Ambient	N/A
Samples	Retention	10 g	4	N/A
	Quality Control	10 g	1	
	Microbiology	10 g	1	