



QUALITY CONTROL CERTIFICATE

Cat. No. MAB194P

Monoclonal Antibody to Human Osteopontin, N-terminal

Description Human Osteopontin (OPN) is a negatively charged hydrophilic protein of 314 amino acids and is subject to significant post translational modifications (PTM) including phosphorylation and glycosylation. Due to its acidic nature and PTMs, OPN runs anomalously by SDS-PAGE. Although its mass is 35kD, apparent molecular weights may range up to 75kD. In addition, OPN is subject to proteolytic modification into smaller molecular weight fragments.

Clone Determination 2F10

Host Mouse

Immunogen Recombinant full-length human Osteopontin

Immunoglobulin Class IgG1 κ

Specificity Specific for the N-terminal fragment of full-length

Osteopontin. Recognizes native human Osteopontin from

purified breast milk and recombinant murine OPN

expressed in E.coli.

Purification Protein A chromatography from ascites

Lot Number

% Purity > % pure determined by size-exclusion chromatography

pl 7.1-7.8

Concentration

Buffer 0.15M PBS

Preservatives None

Expiration date

Quantity

Volume

Applications ELISA, WB

Working dilution Determined by user

Storage Store at -20°C to -80°C; avoid repeat freeze- thaw cycles

References: Alicia Plumer, Hongyi Duan, Sripriya Subramaniam, F. Lee Lucas, Susan Meisfeldt, Ah-Kau Ng, Lucy Liaw. Development of Fragment Specific Osteopontin Antibodies and ELISA for Quantification in Human Metastatic Breast Cancer. BMC Cancer 2008, 8:38.

FOR RESEARCH AND IN VITRO USE ONLY

MAINE BIOTECHNOLOGY SERVICES, INC.