

## Gold Conjugate –Data Sheet

<b>Conjugate:</b>	Goat anti Human IgG (H+L):	40nm Gold
<b>Product Code:</b>	BA.GAHL40s	
<b>Lot Number:</b>	25050089	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
<b>Optical Density at 520nm at manufacture:</b>	QCP064	3.0 – 3.4	3.1	Pass
<b>Mean Particle Diameter:</b>	SPEC-007	37.0 – 43.0nm	40.5nm	Pass
<b>Coefficient of Variation of Particle Size:</b>	SPEC-007	≤ 8.0%	≤ 8.0%	Pass
<b>Buffer:</b>	Conjugate is suspended in 2mM sodium Tetraborate at pH 7.2 containing 0.095% Sodium Azide.			
<b>Stability of this product is guaranteed when stored at 2-8°C until:</b>	22 AUG 2026			

Ref: DS043 Issue 7  
 Print Date: 30 Apr. 25  
 Effective Date: Refer to QPulse

Page 1 of 1

USA: 1-207-517-0284 | International: +44 (0)1495 363 000 | [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)