

## Product Analysis

### EM Goat anti-Rabbit IgG: 15nm Gold

## Product Information

Product Code:	EM.GAR15
Batch Number:	*****
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	***µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 5mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

## Quality Control

Protein Detected:	<***pg	<b>Passed</b>
Optical Density:	*****	<b>Passed</b>
Mean Diameter:	*nm	<b>Passed</b>
Coefficient of Variation:	***%	<b>Passed</b>
Percentage of Single Particles:	**%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**