

CERTIFICATE OF ANALYSIS

POLYCLONAL ANTIBODIES

Ferritin Polyclonal Antibody

| | |
|-----------------|--------------------------|
| Code No. | BP230-3 |
| Grade | Affinity Purified |
| Lot No. | |

RECEIVER INFORMATION

Expiry Date

Storage Temperature 2-8°C recommended, but can be stored frozen (see storage notes below)

Storage Notes **Do not subject to freeze thaw cycles.** If the product is frozen for long term storage this should be done by dispensing into single use aliquots before freezing, also avoiding frost free freezers, which may subject the product to multiple freeze-thaw cycles.

Shipping Temperature 2-8°C

PRODUCT INFORMATION

Immunogen Human Liver Ferritin.

Specificity Reacts with human ferritin

Host Rabbit

Purification Immunoaffinity purification against the antigen.

Presentation Purified polyclonal antibody (Pab) supplied liquid in a PBS buffer, pH7.2 containing 0.09% NaN₃

Applications IP, ELISA capture, Blotting.

HEALTH AND SAFETY

Application For in-vitro use and research or further processing only.

Precaution Potentially hazardous. For specific precautions, see Safety Data Sheet.

Material Safety For further information and technical details, please download a Safety Data Sheet at www.bbisolutions.com or contact your BBI Account Manager.

ANALYSIS

| TESTS | SPECIFICATIONS | RESULTS |
|-------|----------------|---------|
|-------|----------------|---------|

Determination Method:

| | | |
|------------------------------|---------------|-------|
| Protein Concentration | Report result | mg/ml |
|------------------------------|---------------|-------|

| | | | |
|----------------|--|------------------|--|
| Name: | | Position: | |
| Signed: | | Date: | |