

CERTIFICATE OF ANALYSIS

POLYCLONAL ANTIBODIES

Goat anti-Biotin Polyclonal Antibody

Code No.	BP436-3
Grade	Affinity Purified
Lot No.	

RECEIVER INFORMATION

Manufacture Date

Storage Temperature 2-8°C

Shipping Temperature 2-8°C

PRODUCT INFORMATION

Specificity Reacts with Biotin

Host Goat

Purification Affinity purified.

Presentation Purified polyclonal antibody (Pab), supplied lyophilized. Reconstitute by adding 1ml of water. The resulting solution will have the following composition: 10mM phosphate, pH 7.8, 0.15M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Application(s) This product is recommended for use in Western blots and immunohistochemistry (IHC).

HEALTH AND SAFETY

Application For in-vitro use and research or further processing only.

Precaution Potentially hazardous. For specific precautions, see Safety Data Sheet.

Material Safety For further information and technical details, please download a Safety Data Sheet at www.bbisolutions.com or contact your BBI Account Manager.

ANALYSIS

TESTS	SPECIFICATIONS	RESULTS
-------	----------------	---------

Determination Method

Dispensed liquid on concentration value and lyophilised	≥ 85%	Pass
---	-------	------

Name:		Position:	
Signed:		Date:	