

QUALITY CONTROL CERTIFICATE

Cat. No. MAB321P
Dengue NS1 Monoclonal Antibody

Host	Mouse
Clone Determination	10E01
Immunogen	DENV 2 NS1 Recombinant Protein
Immunoglobulin Class	IgG1-kappa
Specificity	Dengue NS1: serovar 1, 2, & 3 recognition
Purification	Protein A chromatography from <i>in vitro</i> production
Lot Number	
% Purity	> 95% pure determined by size-exclusion chromatography
pI	6.0 – 7.0
Concentration	
Buffer	1 x PBS pH 7.6
Preservatives	None
Expiration Date	5 years from date of manufacture
Quantity	
Volume	
Applications	ELISA pairs with MAB320P, suitable for Lateral Flow
Storage	-20°C

FOR RESEARCH AND *IN VITRO* USE ONLY

MAINE BIOTECHNOLOGY SERVICES, INC.
1037R FOREST AVENUE, PORTLAND, ME 04103 USA ~ Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595
Email: sales@mainebiotechnology.com ~ www.mainebiotechnology.com