

## Product Analysis

### EM Rabbit anti Goat IgG (H&L): 10nm Gold

#### Product Information

Product Code:	EM.RAG10
Batch Number:	*****
Protein Absorption:	Not Absorbed
Protein Concentration:	***µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

#### Quality Control

Protein Detected:	<***pg after silver enhancement	<b>Passed</b>
Optical Density:	*****	<b>Passed</b>
Mean Diameter:	*nm	<b>Passed</b>
Coefficient of Variation:	***%	<b>Passed</b>
Percentage of Single Particles:	**%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**