

HD Grade Data Sheet - 20nm Gold Colloid

Information	
Product	Gold colloid – 20nm
Product Code	HD.GC20.OD5
Batch Number	XXXXXXX
Storage Buffer	Suspended in H2O, no preservative Residual chemical left from manufacture
Storage Instructions	Store at 2-8°C – DO NOT FREEZE
Date of Manufacture	Date Month Year (Use within one month of opening)

Technical Data	
Mean Diameter	XX.Xnm
Coefficient of Variation	≤8%
Optical Density @ 520nm	XX.XX
Gold Chloride Concentration	0.05%

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483