

RNP/Sm (non recombinant; bovine)

Antigen Specification

Product Number: 11600

Description:

Bovine RNP/Sm ribonucleoprotein complex.
Natural antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human auto-antibodies.

Origin:

Native. Purified from bovine tissues by protein-chemical methods.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with anti-RNP/Sm autoantibody-positive sample.

Immunological tests/Functionality:

Immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided.

Copyright 2013 - 2021 DIARECT GmbH