

## Ro/SS-A (60 kDa; non recombinant; bovine)

### Antigen Specification

**Product Number: 15500**

**Description:**

Bovine Ro/SS-A protein (60 kDa).  
Natural antigen for in vitro research  
and manufacturing use only.

**Immunological function:**

Binds IgG-type human auto-antibodies.

**Origin:**

Native. Purified from bovine tissues by  
protein-chemical methods.

**Biochemical tests:**

SDS-PAGE (purity > 60%); Western  
blot with antiRo/SS-A 60 kDa  
autoantibody-positive sample.

**Immunological tests/Functionality:**

Standard ELISA test (checkerboard  
analysis of positive/negative samples);  
immunodot analyses with  
positive/negative samples.

**Recommended buffer/storage and  
handling conditions:**

Recommendations for storage buffer:  
neutral to slightly alkaline pH. Storage  
conditions: -70°C or below.  
Repeated freeze/thaw cycles should  
be avoided.

**Coating concentration:**

0.2-0.6 µg/mL (depending on the type  
of ELISA plate and coating buffer).  
Suitable for labeling of functional  
groups.

Copyright 2006 - 2021 DIARECT GmbH