

BPI (non recombinant)

Antigen Specification

Product Number: 19200

Description:

Human BPI (bactericidal permeability-increasing protein).

Natural antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human auto-antibodies.

Origin:

Native. Isolated from human buffy coats.

Biochemical tests:

SDS-PAGE (purity > 80%).

Immunological tests/Functionality:

Immunodot analyses.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: acidic pH and 10% glycerol as cryoprotective agent.

Storage conditions: -70°C or below.

Repeated freeze/thaw cycles should be avoided.

Copyright 2009 - 2021 DIARECT GmbH