

# GP2

## **Antigen Specification**

## Product Number: 19600

#### **Description:**

Human GP2. Recombinant antigen for in vitro research and manufacturing use only.

#### Immunological function:

Binds IgM- and IgG-type human autoantibodies.

#### Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

#### **Expression construct:**

cDNA coding for the free form of the human pancreatic secretory granule membrane major glycoprotein/GP2 protein fused to a deca-histidine purification tag.

#### **Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with monoclonal anti-GP2 antibody.

### Calculated molecular weight:

55 kDa (protein component excluding glycosylation; observed molecular weight approx. 68 kDa).

# Calculated isoelectric point: pH 5.4

### Immunological tests/Functionality:

Standard ELISA test (checkerboard analysis of positive/negative samples).

# Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

#### **Coating concentration:**

0.5-0.9 μg/mL (depending on the type of ELISA plate and coating buffer). Suitable for labeling of functional groups.

Copyright 2010 - 2021 DIARECT GmbH