

## ***Borrelia spielmanii* DbpA**

### **Antigen Specification**

**Product Number: 42900**

**Description:**

*Borrelia spielmanii* (strain A14S) decorin-binding protein A (DbpA). Recombinant antigen for in vitro research and manufacturing use only.

**Immunological function:**

Binds IgG- and IgM-type human antibodies.

**Origin:**

Recombinant. Expressed in *E. coli* bacterial cells.

**Expression construct:**

Full-length cDNA coding for *Borrelia spielmanii* DbpA fused to a deca-histidine purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with patient sample.

**Calculated molecular weight:**

17 kDa

**Calculated isoelectric point:**

pH 6.89

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

**Immunological tests/Functionality:**

Immunodot analyses with Lyme disease positive/negative samples.

**Recommended buffer/storage and handling conditions:**

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2017 - 2021 DIARECT GmbH