

## Ara h 2.0201

### Allergen Specification

**Product Number: 50100**

**Description:**

*Arachis hypogaea* Ara h 2.0201.  
Recombinant allergen for in vitro  
research and manufacturing use only.

**Immunological function:**

Binds IgE-type human antibodies.

**Origin:**

Recombinant. Expressed by  
recombinant baculovirus (*Autographa  
californica* multiple nuclear  
polyhedrosis virus; AcMNPV) infection  
of *Spodoptera frugiperda* Sf9 insect  
cells.

**Expression construct:**

Full-length cDNA coding for Ara h  
2.0201 fused to a deca-histidine  
purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western  
blot with patient sample.

**Calculated molecular weight:**

19 kDa

**Calculated isoelectric point:**

pH 6.47

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte  
Lysate) chromogenic endotoxin assay.

**Immunological tests/Functionality:**

Immunodot analyses with  
positive/negative samples.

**Recommended buffer/storage and  
handling conditions:**

Recommendations for storage buffer:  
neutral to slightly alkaline pH and 6 M  
urea. Storage conditions: -70°C or  
below.

Repeated freeze/thaw cycles should  
be avoided.

Copyright 2016 - 2021 DIARECT GmbH