

## Bet v 1.0101

### Allergen Specification

**Product Number: 50200**

**Description:**

*Betula verrucosa* (*Betula pendula*) Bet v 1.0101.

Recombinant allergen for in vitro research and manufacturing use only.

**Immunological function:**

Binds IgE-type human antibodies.

**Origin:**

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

**Expression construct:**

Full-length cDNA coding for Bet v 1.0101 fused to a deca-histidine purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with patient sample.

**Calculated molecular weight:**

19 kDa

**Calculated isoelectric point:**

pH 6.25

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

**Immunological tests/Functionality:**

Immunodot analyses with positive/negative samples.

**Recommended buffer/storage and handling conditions:**

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2015 - 2021 DIARECT GmbH