

| | |
|--|--|
| Product: | INTRINSIC FACTOR |
| Product code: | IF1 |
| CAS number: | 9008-12-2 |
| EINECS number: | 232-711-9 |
| Alternative name: | Vitamin B12-binding protein, Factor VII, Intrinsic Factor protein. |
| Source: | Porcine stomach pyloric section |
| Form: | A solution in RO water. |
| IF Content: (IF functional ELISA) | > 1.00 x10⁶ ng/ml |
| Purity: (Size exclusion chromatography) | > 95% |
| SDS-Page: | One major band at 55kDa |
| Protein concentration: (Determined by A ₂₈₀ /1.4) | 0.4 – 0.6 mg protein/ml |
| Storage conditions: | Store between -25°C and -10°C. |