

Jo-1 humAb IgG

Antibody Specification

Product Number: 37300

Host: GammaPrim® mouse Clonality: monoclonal Isotype: IgG 1 / kappa

Description:

Human chimeric monoclonal antibody raised against recombinant Jo-1. For in vitro research and manufacturing use only.

Immunogen/Epitope:

Human Jo-1

Immunological tests/Functionality:

ELISA and membrane-based immunoassay with DIARECT's Jo-1 antigen.

Origin:

Produced by heterohybridoma cell line (constructed by fusing peripheral B lymphocytes from a GammaPrim® mouse immunized with Jo-1 to a hetero-hybridoma fusion cell line).

Concentration for assay use:

Concentrations for use as standard and/or positive control reagent in immunodiagnostic assays have to be individually determined by user.

Recommended buffer, storage and handling conditions:

Recommendations for storage buffer: neutral pH and 10% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided.

Copyright 2020 - 2021 DIARECT GmbH