

Product Analysis

EM Goat anti-Rabbit IgG: 40nm Gold

Product Information

Product Code:	EM.GAR40
Batch Number:	010056
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	4.22µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 4°C

Quality Control

Protein Detected:	<400pg	Passed
Optical Density:	5.0	Passed
Mean Diameter:	41.4nm	Passed
Coefficient of Variation:	5.6%	Passed
Percentage of Single Particles:	95%	Passed

For in vitro research only. Not for diagnostic use