

Product Analysis

EM Goat anti-Guinea Pig IgG: 10nm Gold

Product Information

Product Code:	EM.GAG10
Batch Number:	16498
Protein Absorption:	Not Absorbed
Protein Concentration:	17.8µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<4pg after silver enhancement	Passed
Optical Density:	3.1	Passed
Mean Diameter:	9.3nm	Passed
Coefficient of Variation:	7.4%	Passed
Percentage of Single Particles:	98%	Passed

For in vitro research only. Not for diagnostic use