

## **Product Analysis**

### **EM Goat anti-Mouse IgG (H&L): 5nm Gold**

## **Product Information**

Product Code:	EM.GMHL5
Batch Number:	16717
Protein Absorption:	Minimum cross reaction to Human, Bovine and Horse serum proteins
Protein Concentration:	15.68µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

## **Quality Control**

Protein Detected:	<4pg after silver enhancement	<b>Passed</b>
Optical Density:	3.1	<b>Passed</b>
Mean Diameter:	5.5nm	<b>Passed</b>
Coefficient of Variation:	10.1%	<b>Passed</b>
Percentage of Single Particles:	94%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.  
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483