

The Immunogold Specialist

Product Analysis

EM Goat anti-Human IgG: 10nm Gold

Product Information

Product Code: EM.GAH10
Batch Number: 10071
Protein Absorption: Not absorbed
Protein Concentration: 18.7 µg / ml
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA
15mM NaN₃ (preservative) and 30% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected:	7pg after silver enhancement	Passed
Optical Density:	3.0	Passed
Mean Diameter:	10.4nm	Passed
Coefficient of Variation:	7.1%	Passed
Percentage Single Particles:	96%	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use