

Gold Conjugate –Data Sheet

Conjugate:	Goat anti Human IgG:	40nm Gold
Product Code:	BA.GAH40	
Lot Number:	25060117	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
Optical Density at 520nm at manufacture:	QCP064	10.0 – 10.4	10.3	Pass
Mean Particle Diameter:	SPEC-007	37.0 – 43.0nm	40.6nm	Pass
Coefficient of Variation of Particle Size:	SPEC-007	≤ 8.0%	≤ 8.0%	Pass
Buffer:	Conjugate is suspended in 2mM sodium Tetraborate at pH 8.2 containing 0.095% Sodium Azide.			
Stability of this product is guaranteed when stored at 2-8°C until:	09 Jun 2027			

Ref: DS013 Issue 10
 Print Date: 09 Jun. 25
 Effective Date: Refer to QPulse

Page 1 of 1

USA: 1-207-517-0284 | International: +44 (0)1495 363 000| info@bbisolutions.com | www.bbisolutions.com