

Product Analysis

EM Protein A: 10nm Gold

Product Information

Product Code:	EM.PAG10
Batch Number:	22120079
Protein Absorption:	Not Applicable
Protein Concentration:	5.02µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 30% Glycerol
Storage:	Store at 2-8°C

Quality Control

Protein Detected:	<100pg after silver enhancement	Passed
Optical Density:	3.2	Passed
Mean Diameter:	9.7nm	Passed
Coefficient of Variation:	7.9%	Passed
Percentage of Single Particles:	98%	Passed

For in vitro research only. Not for diagnostic use

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | info@bbisolutions.com | www.bbisolutions.com