

## **Product Analysis**

### **EM Goat anti-Rabbit IgG: 10nm Gold**

## **Product Information**

Product Code:	EM.GAR10
Batch Number:	25070079
Protein Absorption:	Minimum cross reaction to Human serum proteins
Protein Concentration:	13.18µg/ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN <sub>3</sub> (preservative) and 30% Glycerol
Storage:	Short Term Storage at 2-8°C (<6months) Long Term Storage -20°C (>6months)

## **Quality Control**

Protein Detected:	<7pg after silver enhancement	<b>Passed</b>
Optical Density:	3.07	<b>Passed</b>
Mean Diameter:	9.6nm	<b>Passed</b>
Coefficient of Variation:	10%	<b>Passed</b>
Percentage of Single Particles:	100%	<b>Passed</b>

**For in vitro research only. Not for diagnostic use**

USA: 1-800-423-8199 | International: +44 (0)2920 767 499 | [info@bbisolutions.com](mailto:info@bbisolutions.com) | [www.bbisolutions.com](http://www.bbisolutions.com)

Part of **BBI Group**

Registered office: c/o Berry Smith LLP, Haywood House, Dumfries Place, Cardiff, CF10 3GA.  
BBI Solutions is the trading name of BBI Solutions OEM Limited. Registered in England and Wales Number 8368483